**Primerdesign**<sup>™</sup> Ltd

# **Aspergillus**

Selective screening kit genesig® kit

150 tests

GENESIG

Kits by Primerdesign

For general laboratory and research use only

## Introduction to Aspergillus

This kit detects four species of Aspergillus; aspergillus fumigatus, aspergillus niger, aspergillus terreus and aspergillus flavus.

Aspergillus fumigatus is a pathogenic fungus of the genus Aspergillus. It is a saprotroph that is found widely in nature, typically in soil and decaying organic matter, such as compost heaps, and is the causative agent for the disease aspergillosis. A. fumigatus has a haploid 29.4 Mb genome that forms eight chromosomes. It forms colonies that produce hyphae known as conidiophores that produce thousands of minute grey-green conidia (2–3 µm) that readily become airborne. Inhaling the spores can result in allergic bronchopulmonary aspergillosis, pulmonary aspergilloma or invasive aspergillosis. Infection occurs after inhalation and it is estimated that everybody inhales several hundred spores each day; typically, these are quickly eliminated by the immune system in healthy individuals. The fungus is more likely to become pathogenic in immunocompromised individuals.

Aspergillus terreus is a filamentous ascomycete fungus with a global distribution, commonly found in soil, compost and dust. It is less common that other species of the Aspergillus genus (e.g. A.fumigatus) however is often the cause of opportunistic aspergillosis. A.terreus was originally thought to be asexual however it is now known to be capable of sexual reproduction. The A.terreus genome is approximately 35 Mbp in size and has a 50-60% GC content. A.terreus has roughly 10,000 protein-coding genes, most of which contain multiple small introns. A.terreus produces characteristic Aspergillus aerial hyphae as well as lateral cells (aleurospores) without conidiophore structures in submerged culture. Invasive aspergillosis as a result of A.terreus infection causes significantly higher mortality than infection with other Aspergillus species.

Aspergillus niger is one of the most common species of the Aspergillus genus. A.niger is a common contaminant of food causing a disease called black mould on certain fruits and vegetables including grapes, apricots, peanuts and onions. As with other species of Aspergillus, A.niger is ubiquitous in the soil. A.niger can cause aspergillosis however infections in humans are generally very rare.

Aspergillus flavus is a common contaminant of cereal grains, tree nuts and legumes. Infections can occur pre-harvest however often the plant shows no symptoms until post-harvest storage and/or during transport. Many strains produce significant quantities of toxic compounds known as mycotoxins, which are toxic to mammals when consumed. Like many other species of Aspergillus, A. flavus can cause Aspergillosis in immunocompromised individuals.

# **Specificity**

The kit contains individual primers directed against the four following *aspergillus* species; *aspergillus fumigatus*, *aspergillus niger*, *aspergillus terreus* and *aspergillus flavus*.

The primers and probe sequences in this kit have 100% homology with a broad range of clinically relevant reference sequences based on a comprehensive bioinformatics analysis. Please note, this assay may also detect strains of *aspergillus oryzae*.

If you require further information or have a specific question about the detection profile of this kit, then please send an email to <a href="mailto:enquiry@primerdesign.co.uk">enquiry@primerdesign.co.uk</a> and our bioinformatics team will answer your question.

### Kit contents

Aspergillus\_SCRN primer/probe mix (150 reactions BROWN)
 FAM and VIC labelled (see table below)

Target	Fluorophore
Aspergillus fumigatus	FAM
Aspergillus niger	FAM
Aspergillus terreus	FAM
Aspergillus flavus	FAM
Internal extraction control	VIC

- Aspergillus\_SCRN positive control template (RED)
   Positive signal for FAM channel
- Internal extraction control DNA (150 reactions BLUE)
- RNase/DNase free water (WHITE) for resuspension of primer/probe mix
- Template preparation buffer (YELLOW)
   for resuspension of positive control template and internal control template

### Reagents and equipment to be supplied by the user

#### **Real-Time PCR Instrument**

Ability to detect FAM and VIC channels

#### **DNA** extraction kit

This kit is recommended for use with genesig Easy DNA/RNA extraction kit. However, it is designed to work well with all processes that yield high quality DNA with minimal PCR inhibitors.

### oasig<sup>™</sup> lyophilised or Precision<sup>®</sup>PLUS 2X qPCR Master Mix

This kit is intended for use with oasig or PrecisionPLUS 2X qPCR Master Mix.

### **Pipettors and Tips**

**Vortex and Centrifuge** 

Thin walled 1.5ml PCR reaction tubes

## Kit storage and stability

This kit is stable at room temperature but should be stored at -20°C on arrival. Once the lyophilised components have been resuspended they should not be exposed to temperatures above -20°C for longer than 30 minutes at a time and unnecessary repeated freeze/thawing should be avoided. The kit is stable for six months from the date of resuspension under these circumstances.

Primerdesign does not recommend using the kit after the expiry date stated on the pack.

# Suitable sample material

All kinds of sample material suited for PCR amplification can be used. Please ensure the samples are suitable in terms of purity, concentration, and RNA/DNA integrity. Always run at least one negative control with the samples. To prepare a negative control, replace the template DNA sample with RNase/DNase free water.

# **Dynamic range of test**

Under optimal PCR conditions genesig kits have very high priming efficiencies of >90% and can detect between 1x10<sup>8</sup> and 1x10<sup>2</sup> copies of target template.

### **Notices and disclaimers**

This product is developed, designed and sold for research purposes only. It is not intended for human diagnostic or drug purposes or to be administered to humans unless clearly expressed for that purpose by the Food and Drug Administration in the USA or the appropriate regulatory authorities in the country of use. During the warranty period Primerdesign genesig detection kits allow precise and reproducible data recovery combined with excellent sensitivity. For data obtained by violation to the general GLP guidelines and the manufacturer's recommendations the right to claim under guarantee is expired. PCR is a proprietary technology covered by several US and foreign patents. These patents are owned by Roche Molecular Systems Inc. and have been sub-licensed by PE Corporation in certain fields. Depending on your specific application you may need a license from Roche or PE to practice PCR. Additional information on purchasing licenses to practice the PCR process may be obtained by contacting the Director of Licensing at Roche Molecular Systems, 1145 Atlantic Avenue, Alameda, CA 94501 or Applied Biosystems business group of the Applera Corporation, 850 Lincoln Centre Drive, Foster City, CA 94404. In addition, the 5' nuclease assay and other homogeneous amplification methods used in connection with the PCR process may be covered by U.S. Patents 5,210,015 and 5,487,972, owned by Roche Molecular Systems, Inc, and by U.S. Patent 5,538,848, owned by The Perkin-Elmer Corporation.

### **Trademarks**

Primerdesign™ is a trademark of Primerdesign Ltd.

genesig® is a registered trademark of Primerdesign Ltd.

oasig  $\ensuremath{^{\text{TM}}}$  is a trademark of Primerdesign Ltd.

The PCR process is covered by US Patents 4,683,195, and 4,683,202 and foreign equivalents owned by Hoffmann-La Roche AG. BI, ABI PRISM® GeneAmp® and MicroAmp® are registered trademarks of the Applera Genomics (Applied Biosystems Corporation). BIOMEK® is a registered trademark of Beckman Instruments, Inc.; iCycler™ is a registered trademark of Bio-Rad Laboratories, Rotor-Gene is a trademark of Corbett Research. LightCycler™ is a registered trademark of the Idaho Technology Inc. GeneAmp®, TaqMan® and AmpliTaqGold® are registered trademarks of Roche Molecular Systems, Inc. The purchase of the Primerdesign reagents cannot be construed as an authorization or implicit license to practice PCR under any patents held by Hoffmann-LaRoche Inc.

## Principles of the test

#### **Real-time PCR**

Individual primers and probes designed for each aspergillus species have been combined into a single reaction and these can be detected through the FAM channel as described in the kit contents.

The primer/probe mix provided exploits the so-called TagMan® principle. During PCR amplification, forward and reverse primers hybridise to the target DNA. Fluorogenic probes are included in the same reaction mixture, which consists of a DNA probe labelled with a 5'-dye and a 3'-quencher. During qPCR amplification, the probe is cleaved, and the reporter dye and quencher are separated. The resulting increase in fluorescence can be detected on a range of qPCR platforms.

#### Positive control

For a positive control, the kit contains a single positive control that contains template for one of the aspergillus species and will give a result through the FAM channel. Each time the kit is used, at least one positive control reaction must be included in the run. A positive result indicates that the primer/probe mix and master mix are working properly in that run. If a negative result is obtained the test results are invalid and must be repeated. Care should be taken to ensure that the positive controls do not contaminate any other kit component, which would lead to false positive results. This can be achieved by handling these components in a post PCR environment. Care should also be taken to avoid cross-contamination of other samples when adding the positive control to the run. This can be avoided by sealing all other samples and negative controls before pipetting the positive control into the positive control well.

### **Negative control**

To confirm the absence of contamination, a negative control reaction should be included every time the kit is used. For this reaction, the RNase/DNase free water should be used instead of template. A negative result indicates that the reagents have not become contaminated while setting up the run.

#### **Internal DNA extraction control**

When performing DNA extraction, it is often advantageous to have an exogenous source of DNA template that is spiked into the lysis buffer. This control DNA is then co-purified with the sample DNA and can be detected as a positive control for the extraction process. Successful co-purification and qPCR for the control DNA also indicates that PCR inhibitors are not present at a high concentration. The Aspergillus\_SCRN primer/probe mix supplied already contains separate primers and probe to detect the exogenous DNA using qPCR. The primers are present at PCR limiting concentrations which allows multiplexing with the target sequence primers. Amplification of the control DNA does not interfere with detection of the Aspergillus target DNA even when present at low copy number. The Internal control is detected through the VIC channel and gives a Cq value of 28+/-3.

# **Resuspension protocol**

To minimise the risk of contamination with foreign DNA, we recommend that all pipetting be performed in a PCR clean environment. Ideally this would be a designated PCR lab or PCR cabinet. Filter tips are recommended for all pipetting steps.

1. Pulse-spin each tube in a centrifuge before opening.

This will ensure lyophilised primer/probe mix is in the base of the tube and is not spilt upon opening the tube.

2. Resuspend the Aspergillus\_SCRN primer/probe mix in the RNase/DNase free water supplied, according to the table below:

To ensure complete resuspension, vortex the tube thoroughly.

Component – resuspend in water	
Pre-PCR pack	
Aspergillus_SCRN primer/probe mix (BROWN)	165 µl

3. Resuspend internal control template and positive control template and in the template resuspension buffer supplied, according to the table below:

To ensure complete resuspension, vortex the tube thoroughly.

Component – resuspend in template preparation buffer	
Pre-PCR heat-sealed foil	
Internal extraction control DNA (BLUE)	600 µl
Post-PCR heat-sealed foil	
Positive control template (RED)*	500 µl

\* This component contains high copy number template and is a VERY significant contamination risk. It must be opened and handled in a separate laboratory environment, away from the other components.

### **DNA** extraction

The internal extraction control DNA can be added either to the DNA lysis/extraction buffer or to the DNA sample once it has been resuspended in lysis buffer.

DO NOT add the internal extraction control DNA directly to the unprocessed biological sample as this will lead to degradation and a loss in signal.

- 1. Add 4µl of the Internal extraction control DNA (**BLUE**) to each sample in DNA lysis/extraction buffer per sample.
- 2. Complete DNA extraction according to the manufacturer's protocols.

# qPCR detection protocol

1. For each DNA sample prepare a reaction mix according to the table below: Include sufficient reactions for positive and negative controls.

Component	Volume
oasig or PrecisionPLUS 2X qPCR Master Mix	10 µl
Aspergillus_SCRN primer/probe mix (BROWN)	1 µl
RNase/DNase free water (WHITE)	4 µl
Final volume	15 µl

- 2. Pipette 15µl of these mixes into each well according to your qPCR experimental plate set up.
- 3. Prepare sample DNA templates for each of your samples.
- 4. Pipette 5µl of DNA template into each well according to your experimental plate set up.

For negative control (NTC) wells use 5µl of RNase/DNase free water. The final volume in each well is 20µl.

5. Pipette 5µl of positive control template into the positive control wells.

## qPCR amplification protocol

Amplification conditions using oasig or PrecisionPLUS 2X qPCR Master Mix

	Step	Time	Temp
	Enzyme activation	2 mins	95°C
X 50 cycles	Denaturation	10 secs	95°C
	DATA COLLECTION*	60 secs	60°C


<sup>\*</sup> Fluorogenic data should be collected during this step through the FAM and VIC channels.

# Interpretation of results


Target (FAM)	Internal control (VIC)	Positive Control	Negative Control	Interpretation
+	+/-	+	-	POSITIVE QUALITATIVE RESULT
-	+	+	-	NEGATIVE RESULT
+/-	+/-	+	≤35	EXPERIMENT FAILED  Due to test contamination
+/-	+/-	+	>35	*
-	-	+	-	SAMPLE PREPARATION FAILED
+/-	+/-	-	+/-	EXPERIMENT FAILED

Positive control template (**RED**) is expected to amplify between Cq 16 and 23. Failure to satisfy these quality control criteria is a strong indication that the experiment has been compromised.

\* Where the test sample is positive, and the negative control is positive with a Cq >35, the sample must be reinterpreted based on the relative signal strength of the two results:



If the sample amplifies > 5 Cq earlier than the negative control, then the sample should be reinterpreted (via the table above) with the negative control verified as negative.



If the sample amplifies < 5 Cq earlier than the negative control, then the positive sample result is invalidated, and the result should be determined inconclusive due to test contamination. The test for this sample should be repeated.

### **Internal PCR control**

The Cq value obtained with the internal control will vary significantly depending on the extraction efficiency, the quantity of DNA added to the PCR reaction and the individual machine settings. Cq values of 28±3 are within the normal range. When amplifying an Aspergillus\_SCRN sample with a high genome copy number, the internal extraction control may not produce an amplification plot. This does not invalidate the test and should be interpreted as a positive experimental result.

#### **Positive control**

The positive control well should give an amplification plot through the FAM channel. The positive control signals indicate that the kit is working correctly. If a negative result is obtained in the FAM channel, then the test results are invalid and must be repeated

#### No template control (NTC)

The NTC should give a flat line (flat amplification plots). Signals in the NTC indicate cross contamination during plate set up.

### Sample data

Presence of Aspergillus species is detected in the FAM channel as indicated in the kit contents section. A positive signal in the FAM channel therefore indicates the presence of any of the four Aspergillus species detected by this kit, this kit does not allow you to determine which specific species are present in your samples. Mixed infections, where more than one species are present in the same sample, may occur. In these situations, detection of mixed infections can lead to unusually shaped curves however this is normal, and the data should not be interpreted any differently.